
KNOW MEDICAL

Kit Anestesia

Kits for Anaesthesia

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Knowmedical has introduced between own products the epidural and spinal anaesthesia.

The epidural anaesthesia is an anaesthesia of rachis through an injection of local anaesthetic drug in the extradural space.

The introduction in the epidural space could be made in the cervical, thoracic, lumbar and sacral area.

The epidural space

This area extends from the magnum foramen, where arrive the dura pleura, to the sacra coccyx ligament, which includes the sacral iatum. This is an open space, because is crossed from 31 pair of spinal nerves and communicates with the Para vertebral spaces.

The epidural space in lumbar region, in the introduction point of the needle on the median line, is deep about mm. 4-6.

The distance of the epidural space from the skin changes from cm. 2,5 to cm. 9 cm with medium value of cm. 4-5.

Tuohy needle is the most used for the puncture in the epidural space. The size changes between 16G to 20G. The edge is curved towards the height to allow the exit of catheter and determinate his direction, moreover it is less cutting for reducing the risk of dura puncture. The needle has to be graduated in cm. to guarantee the measure of the reached depth.

Each needle has a *mandrin* that run easy in order to avoid the obstruction of lumen of the needle from little pieces of skin and tissue. This element is taken away before the introduction of the yellow ligament to allow the application of technique that identifies the epidural space (Technique of loss of resistance and Technique of pending drop).

The *catheter* is in nylon to allow obtaining a product neither too soft (risk to clinking) nor too rigid (risk of wounds). Moreover the radiopaque line has the function to see it in case of breaking.

The epidural catheter must have these characteristics:

- Not to be either too soft or too rigid.
- Must be resistant for not breaking itself during his extraction, and for not clinking during the puncture.
- The diameter must be enough reduced in order to allow his passage into the Tuohy needle.
- The inner diameter must allow the injection of anaesthetic drug without too resistance.
- The length must be sufficient to fix it easily in the thoracic region.
- The edge must be soft and regular to avoid venous perforation and of the dura mother.
- It must have the cm. to allow the exact identification of the length of catheter which is in the epidural space.
- It could have only one hole, distal or side.

The *syringe* is a device that allows identifying the good position of the needle in the epidural space. The syringe could be two different models:

- L.O.R. syringe (loss of resistance) to low resistance
- Standard syringe with best sensibility

In order to better identify the epidural space it is necessary that the plunger move easily into the syringe body.


The *adaptor* has a male Luer Lock cap that allow both to maintain the sterilization inside until the use even if the package is open, and the correct position of catheter.

The *filter* is specific for the administration of anaesthetics in epidural with the aim to cancel the risks of particles and/or contamination. It has a male Luer Lock cap e a vented cap that allow to the sterilization inside until the use even if the package is open.

It could be standard with porosity 0,2µ, 1,2µ, 5µ or 0,2µ for bacterial retention of endotoxine.

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Epidural Kit for continuous anaesthesia composed by:

- Tuohy needle
- Epidural catheter
- Epidural filter 0,22 μ
- LOR Syringe 10 ml Luer-Lock
- Tuohy borst adapter
- Luer slip Connector

CODICE CODE	TUOHY	CATETERE CATHETER	IMBALLO INTERNO INNER PACKAGING	IMBALLO ESTERNO OUTER PACKAGING
EP0 016 O	16 G	3fr x 88 cm	20	120
EP0 017 O	17 G	3fr x 88 cm	20	120
EP0 018 O	18 G	20 G	20	120

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Also you have the possibility of customize kit with accessories, according to the needs and request of the customers (drapes, gloves, gauzes, klemmer, clamps, etc.)

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

Quincke Spinal needles with introducer

LUNGHEZZA 90mm. LENTGH 90mm.	
CODICE CODE	MISURA SIZE
SP2 216 O	16Gx90mm
SP2 217 O	17Gx90mm
SP2 218 O	18Gx90mm
SP2 219 O	19Gx90mm
SP2 220 O	20Gx90mm
SP2 221 O	21Gx90mm
SP2 222 O	22Gx90mm
SP2 223 O	23Gx90mm
SP2 224 O	24Gx90mm
SP2 225 O	25Gx90mm
SP2 226 O	26Gx90mm
SP2 227 O	27Gx90mm

LUNGHEZZA 103mm. LENTGH 103mm.	
CODICE CODE	MISURA SIZE
SP2 316 O	16Gx103mm
SP2 317 O	17Gx103mm
SP2 318 O	18Gx103mm
SP2 319 O	19Gx103mm
SP2 320 O	20Gx103mm
SP2 321 O	21Gx103mm
SP2 322 O	22Gx103mm
SP2 323 O	23Gx103mm
SP2 324 O	24Gx103mm
SP2 325 O	25Gx103mm
SP2 326 O	26Gx103mm
SP2 327 O	27Gx103mm



LUNGHEZZA 120mm. LENTGH 120mm.	
CODICE CODE	MISURA SIZE
SP2 416 O	16Gx120mm
SP2 417 O	17Gx120mm
SP2 418 O	18Gx120mm
SP2 419 O	19Gx120mm
SP2 420 O	20Gx120mm
SP2 421 O	21Gx120mm
SP2 422 O	22Gx120mm
SP2 423 O	23Gx120mm
SP2 424 O	24Gx120mm
SP2 425 O	25Gx120mm
SP2 426 O	26Gx120mm
SP2 427 O	27Gx120mm

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
Pencil Point Spinal needles with introducer

LUNGHEZZA 90mm. <i>LENTGH 90mm.</i>	
CODICE <i>CODE</i>	MISURA <i>SIZE</i>
SP3 222 O	22Gx90mm
SP3 223 O	23Gx90mm
SP3 224 O	24Gx90mm
SP3 225 O	25Gx90mm
SP3 226 O	26Gx90mm
SP3 227 O	27Gx90mm

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KNOWMEDICAL S.r.l.